

Package: axR (via r-universe)

July 16, 2026

Title Serial Communication and Data Retrieval for Axivity Devices

Version 0.1.2

Description Serial communication and data retrieval for Axivity AX3/AX6 accelerometer devices: device discovery, status, settings, data download, and reading recorded '.cwa'/AX6 binary files. Wraps the Open Movement Project's 'OMAPI' C library (vendored in src/omapi, BSD 2-clause, Newcastle University; see src/omapi/LICENSE.TXT), rather than reimplementing the serial protocol or binary file format directly. On Linux, device discovery requires 'libudev-dev' (Debian/Ubuntu) or 'systemd-devel' (Fedora/RHEL) to be installed separately -- see SystemRequirements. macOS and Windows use OS frameworks/APIs directly, with no separate install step. Part of the Circadia Lab R ecosystem.

License MIT + file LICENSE

Encoding UTF-8

Roxygen list(markdown = TRUE)

Depends R (>= 4.1.0)

Imports Rcpp (>= 1.0.0)

LinkingTo Rcpp

SystemRequirements libudev

Suggests covr, pkgdown, testthat (>= 3.1.4), tibble, knitr, rmarkdown

VignetteBuilder knitr

Config/testthat/edition 3

Lifecycle experimental

URL <https://axr.circadia-lab.uk>, <https://github.com/circadia-bio/axR>

BugReports <https://github.com/circadia-bio/axR/issues>

Config/roxygen2/version 8.0.0

RoxygenNote 7.3.3

Repository <https://circadia-bio.r-universe.dev>

Date/Publication 2026-07-16 12:14:51 UTC
RemoteUrl <https://github.com/circadia-bio/axR>
RemoteRef v0.1.2
RemoteSha 893b154bed8abe50f0a11ce1580df39dc2915606

Contents

axivity_copy_data	2
axivity_discover	3
axivity_download	4
axivity_download_cancel	4
axivity_download_status	5
axivity_enable_debug_log	5
axivity_get_accel_config	6
axivity_get_accelerometer	7
axivity_get_battery	7
axivity_get_data_info	8
axivity_get_delays	8
axivity_get_ecc	9
axivity_get_memory_health	9
axivity_get_metadata	10
axivity_get_session_id	10
axivity_get_time	11
axivity_is_locked	11
axivity_read_cwa	12
axivity_reset	13
axivity_self_test	13
axivity_send_command	14
axivity_set_led	14

Index	15
--------------	-----------

axivity_copy_data	<i>Copy recorded data directly off a mounted Axivity volume</i>
-------------------	---

Description

A fallback for when `axivity_discover()`/OMAPI device access isn't working, but the device's USB mass-storage volume mounts and is visible in Finder/diskutil regardless – a common split, since the storage side and OMAPI's own IOKit-level device discovery are independent paths (see NEWS.md for the discovery issues hit so far). This bypasses OMAPI and device_id entirely: it's a plain file copy, nothing more. axR doesn't parse .cwa contents – see mrpheus or zeitR for that.

Usage

```
axivity_copy_data(
  device_path,
  dest_dir,
  pattern = "\\\\.cwa$",
  overwrite = FALSE
)
```

Arguments

device_path	Character. Path to the mounted device volume, e.g. <code>"/Volumes/CWA17_46171"</code> . Find this in Finder, or with <code>list.files("/Volumes")</code> .
dest_dir	Character. Destination directory for downloaded files. Created (recursively) if it doesn't already exist.
pattern	Character. Regex to filter which files are copied, matched case-insensitively. Default <code>"\\.cwa\$"</code> .
overwrite	Logical. Overwrite existing files at the destination. Default FALSE.

Value

Character vector of destination file paths that were successfully copied.

axivity_discover	<i>Discover connected Axivity devices</i>
------------------	---

Description

Wraps OMAPI's `OmGetDeviceIds()` plus a handful of per-device status calls into a single data frame. Unlike serial-port probing, this uses OMAPI's own device discovery – including its platform-specific finder (IOKit/DiskArbitration on macOS, SetupAPI on Windows, udev on Linux) – so it should behave the same way OmGui does on the same machine.

Usage

```
axivity_discover()
```

Value

A data frame with one row per connected device: `device_id`, `serial`, `port`, `path`, `firmware_version`, `hardware_version`, `battery_level`. Zero rows if no devices are connected.

axivity_download	<i>Download recorded data off an Axivity device</i>
------------------	---

Description

Wraps OMAPI's `OmBeginDownloading()`, which runs the download on a background thread inside the library. `axR` doesn't parse `.cwa` contents – see `mrpheus` or `zeitR` for that.

Usage

```
axivity_download(
  device_id,
  dest_file,
  offset_blocks = 0L,
  length_blocks = -1L,
  blocking = TRUE
)
```

Arguments

<code>device_id</code>	Integer device ID, as returned by <code>axivity_discover()</code> .
<code>dest_file</code>	Character. Destination file path (existing files at this path are truncated).
<code>offset_blocks</code>	Integer. Start offset of the download, in blocks. Default 0.
<code>length_blocks</code>	Integer. Length to download, in blocks. Default -1 (all).
<code>blocking</code>	Logical. If TRUE (default), block until the download completes, fails, or is cancelled – equivalent to calling <code>axivity_download_wait()</code> immediately after. If FALSE, return as soon as the download starts; poll with <code>axivity_download_status()</code> .

Value

If `blocking = TRUE`, a list with status (one of "complete", "error", "cancelled") and value (a diagnostic code if status is "error"). If `blocking = FALSE`, invisibly NULL.

axivity_download_cancel	<i>Cancel an in-progress download from an Axivity device</i>
-------------------------	--

Description

Cancel an in-progress download from an Axivity device

Usage

```
axivity_download_cancel(device_id)
```

Arguments

device_id Integer device ID, as returned by `axivity_discover()`.

axivity_download_status

Check or wait for an Axivity device's download progress

Description

Check or wait for an Axivity device's download progress

Usage

```
axivity_download_status(device_id)
```

```
axivity_download_wait(device_id)
```

Arguments

device_id Integer device ID, as returned by `axivity_discover()`.

Value

A list with status ("none", "error", "progress", "complete", or "cancelled") and value (percentage complete if status is "progress", a diagnostic code if "error").

axivity_enable_debug_log

Enable OMAPI's internal debug log

Description

A diagnostic escape hatch for when `axivity_discover()` isn't finding a device you expect it to. OMAPI logs internally via its own `OmLog()`, but axR's default log target is NULL (so compiled code doesn't write to stderr unprompted – see NEWS.md). This re-enables it.

Usage

```
axivity_enable_debug_log(file = NULL)
```

Arguments

file Character path to a log file, or NULL (default) to log to stderr. A file is more reliable for diagnostic purposes: OMAPI's discovery thread logs from a background pthread, and raw stderr writes from a non-R thread don't always reach the R console/terminal depending on the frontend – a file sidesteps that ambiguity.

Details

Important: this only controls where log lines go. Whether anything is logged *at all* is controlled by OMAPI's debug level, which is read from the OMDEBUG environment variable once, at `OmStartup()` time – i.e. before `library(axR)` runs. If you're not seeing log output after calling this, set OMDEBUG (e.g. `Sys.setenv(OMDEBUG = "3")`) *before* loading axR, in a fresh R session, then call this function and retry.

Value

Invisibly, the OMAPI status code (negative indicates failure, e.g. the file couldn't be opened).

```
axivity_get_accel_config
```

Get or set an Axivity device's accelerometer configuration

Description

Get or set an Axivity device's accelerometer configuration

Usage

```
axivity_get_accel_config(device_id)
```

```
axivity_set_accel_config(device_id, rate, range)
```

Arguments

<code>device_id</code>	Integer device ID, as returned by <code>axivity_discover()</code> .
<code>rate</code>	Sampling rate in Hz (6 = 6.25, 12 = 12.5, 25, 50, 100, 200, 400, 800, 1600, 3200). Negative = low-power mode.
<code>range</code>	Sampling range in +/- G (2, 4, 8, 16).

Value

`axivity_get_accel_config()` returns a list with `rate` (Hz) and `range` (+/- g).

`axivity_get_accelerometer`*Read an Axivity device's current accelerometer values*

Description

Read an Axivity device's current accelerometer values

Usage

```
axivity_get_accelerometer(device_id)
```

Arguments

`device_id` Integer device ID, as returned by `axivity_discover()`.

Value

A named numeric vector `c(x, y, z)` in units of `g` (raw values are in 1/256 `g`, converted here).

`axivity_get_battery` *Query an Axivity device's battery level and health*

Description

Query an Axivity device's battery level and health

Usage

```
axivity_get_battery(device_id)
```

Arguments

`device_id` Integer device ID, as returned by `axivity_discover()`.

Value

A list with `level_pct` (0-99% = charging, 100% = full) and `recharge_cycles` (lower is better).

`axivity_get_data_info` *Get information about an Axivity device's recorded data*

Description

Get information about an Axivity device's recorded data

Usage

```
axivity_get_data_info(device_id)
```

Arguments

`device_id` Integer device ID, as returned by `axivity_discover()`.

Value

A list with `size_bytes`, `filename` (path on the device's own filesystem, not yet downloaded), `block_size`, `offset_blocks`, `num_blocks`, and `start/end` (POSIXct, the time range of the recorded data).

`axivity_get_delays` *Get or set an Axivity device's delayed activation window*

Description

Get or set an Axivity device's delayed activation window

Usage

```
axivity_get_delays(device_id)
```

```
axivity_set_delays(device_id, start, stop)
```

Arguments

`device_id` Integer device ID, as returned by `axivity_discover()`.
`start, stop` Each either a POSIXct (or coercible), `-Inf` (always record from now / OMAPI's zero sentinel), or `Inf` (never record / OMAPI's infinite sentinel).

Value

`axivity_get_delays()` returns a list with `start` and `stop`, each either a POSIXct, `-Inf` (OMAPI's "always"/zero sentinel), or `Inf` (OMAPI's "never"/infinite sentinel).

axivity_get_ecc	<i>Get or set an Axivity device's error-correcting code (ECC) flag</i>
-----------------	--

Description

Get or set an Axivity device's error-correcting code (ECC) flag

Usage

```
axivity_get_ecc(device_id)
```

```
axivity_set_ecc(device_id, enabled)
```

Arguments

device_id	Integer device ID, as returned by axivity_discover() .
enabled	Logical. Enable or disable ECC.

Value

axivity_get_ecc() returns a logical.

axivity_get_memory_health	<i>Query an Axivity device's NAND flash memory health</i>
---------------------------	---

Description

Query an Axivity device's NAND flash memory health

Usage

```
axivity_get_memory_health(device_id)
```

Arguments

device_id	Integer device ID, as returned by axivity_discover() .
-----------	--

Value

A list with spare_blocks (higher is better, 0 = unusable) and status ("ok", "warning", or "error", using OMAP1's documented thresholds).

`axivity_get_metadata` *Get or set an Axivity device's metadata scratch buffer*

Description

Get or set an Axivity device's metadata scratch buffer

Usage

```
axivity_get_metadata(device_id)
```

```
axivity_set_metadata(device_id, metadata)
```

Arguments

<code>device_id</code>	Integer device ID, as returned by axivity_discover() .
<code>metadata</code>	Character. Metadata to store (up to 448 bytes; longer values are truncated by the device). URL-encode first if it needs to preserve non-ASCII characters.

Value

`axivity_get_metadata()` returns a character string, trimmed of trailing padding.

`axivity_get_session_id`
Get or set an Axivity device's session identifier

Description

Get or set an Axivity device's session identifier

Usage

```
axivity_get_session_id(device_id)
```

```
axivity_set_session_id(device_id, session_id)
```

Arguments

<code>device_id</code>	Integer device ID, as returned by axivity_discover() .
<code>session_id</code>	A value to set as the session ID.

Value

`axivity_get_session_id()` returns a numeric (session IDs can exceed R's 32-bit integer range).

axivity_get_time	<i>Get or set an Axivity device's real-time clock</i>
------------------	---

Description

Get or set an Axivity device's real-time clock

Usage

```
axivity_get_time(device_id)
```

```
axivity_set_time(device_id, time)
```

Arguments

device_id	Integer device ID, as returned by axivity_discover() .
time	A POSIXct (or coercible) date/time to set on the device.

Value

axivity_get_time() returns a POSIXct.

axivity_is_locked	<i>Check or set an Axivity device's anti-tamper lock</i>
-------------------	--

Description

Check or set an Axivity device's anti-tamper lock

Usage

```
axivity_is_locked(device_id)
```

```
axivity_set_lock(device_id, code)
```

```
axivity_unlock(device_id, code)
```

Arguments

device_id	Integer device ID, as returned by axivity_discover() .
code	Integer lock code. 0 = no lock; 0xffff is reserved.

Value

axivity_is_locked() returns a list with locked and has_lock_code (logicals).

axivity_read_cwa	<i>Read a .cwa/AX6 binary file into a tibble</i>
------------------	--

Description

Uses OMAPI's own binary file reader (`omapi-reader.c`, vendored from `libomapi`) to parse a .cwa/AX6 recording directly, rather than reimplementing the binary format from scratch. Returns one row per sample – ready to hand to `zeitR`, or any other downstream actigraphy analysis.

Usage

```
axivity_read_cwa(path)
```

Arguments

path Character. Path to a .cwa/AX6 file.

Details

Unlike the rest of `axR`, this function doesn't talk to a live device at all – it works on a file already on disk (e.g. one retrieved with `axivity_copy_data()` or `axivity_download()`), and doesn't require `axivity_discover()` to have found anything.

Value

A tibble (or plain data frame, if the `tibble` package isn't installed) with one row per sample:

timestamp POSIXct, UTC, with sub-second precision

x, y, z Accelerometer readings, in g

gx, gy, gz Gyroscope readings, raw units (only present if the recording has a gyroscope, e.g. AX6 in GA/GAM mode)

mx, my, mz Magnetometer readings, raw units (only present if the recording has a magnetometer, e.g. AX6 in GAM mode)

light Raw light sensor reading

temperature_c Temperature in degrees Celsius. **Unverified, possibly wrong** – OMAPI's conversion (`OM_VALUE_TEMPERATURE_MC`) hardcodes a formula for one specific temperature sensor chip (MCP9700); a comment beside it in the vendored source notes an alternate formula for a different chip (MCP9701), suggesting this may be hardware/revision-specific. Cross-check against OmGui's own reading for the same file before relying on this.

battery_pct Battery percentage, at the time of this block

sample_rate Sampling rate in Hz, at the time of this block

with `device_id`, `session_id`, and metadata attached as attributes. `timestamp`, `x/y/z`, and `device_id` have been verified correct against a real AX3 file (cross-checked `device_id` against `ioreg` and the Axivity config web tool); `temperature_c` has not.

axivity_reset	<i>Erase an Axivity device's data storage and commit settings</i>
---------------	---

Description

Wraps OMAPI's `OmEraseDataAndCommit()`. Staged settings changes (delays, session ID, meta-data, accelerometer config) only take full effect when this is called.

Usage

```
axivity_reset(device_id, level = "quickformat")
```

Arguments

device_id	Integer device ID, as returned by axivity_discover() .
level	One of "none" (commit metadata only, not recommended – can cause a data/metadata mismatch), "delete" (remove and recreate the data file), "quickformat" (recreate the filesystem), or "wipe" (clear all NAND blocks, cleanest but slowest).

axivity_self_test	<i>Run an Axivity device's built-in self-test</i>
-------------------	---

Description

Run an Axivity device's built-in self-test

Usage

```
axivity_self_test(device_id)
```

Arguments

device_id	Integer device ID, as returned by axivity_discover() .
-----------	--

Value

A list with `passed` (logical) and `diagnostic_code` (an opaque, firmware-specific code; 0 means passed).

`axivity_send_command` *Send a raw command to an Axivity device*

Description

An escape hatch wrapping OMAPI's `OmCommand()`, for anything not covered by axR's typed functions. Not generally recommended – OMAPI's own docs note that incorrect use could lead to unspecified behaviour.

Usage

```
axivity_send_command(device_id, command, expected = "", timeout_ms = 2000L)
```

Arguments

<code>device_id</code>	Integer device ID, as returned by <code>axivity_discover()</code> .
<code>command</code>	Character. The command string to send.
<code>expected</code>	Character. The expected response prefix, or "" if not specified.
<code>timeout_ms</code>	Integer. Timeout in milliseconds. Default 2000.

Value

Character. The device's raw response.

`axivity_set_led` *Set an Axivity device's LED colour*

Description

Set an Axivity device's LED colour

Usage

```
axivity_set_led(device_id, colour)
```

Arguments

<code>device_id</code>	Integer device ID, as returned by <code>axivity_discover()</code> .
<code>colour</code>	One of "auto" (default device-controlled behaviour), "off", "blue", "green", "cyan", "red", "magenta", "yellow", "white".

Index

axivity_copy_data, 2
axivity_copy_data(), 12
axivity_discover, 3
axivity_discover(), 2, 4–14
axivity_download, 4
axivity_download(), 12
axivity_download_cancel, 4
axivity_download_status, 5
axivity_download_status(), 4
axivity_download_wait
 (axivity_download_status), 5
axivity_download_wait(), 4
axivity_enable_debug_log, 5
axivity_get_accel_config, 6
axivity_get_accelerometer, 7
axivity_get_battery, 7
axivity_get_data_info, 8
axivity_get_delays, 8
axivity_get_ecc, 9
axivity_get_memory_health, 9
axivity_get_metadata, 10
axivity_get_session_id, 10
axivity_get_time, 11
axivity_is_locked, 11
axivity_read_cwa, 12
axivity_reset, 13
axivity_self_test, 13
axivity_send_command, 14
axivity_set_accel_config
 (axivity_get_accel_config), 6
axivity_set_delays
 (axivity_get_delays), 8
axivity_set_ecc (axivity_get_ecc), 9
axivity_set_led, 14
axivity_set_lock (axivity_is_locked), 11
axivity_set_metadata
 (axivity_get_metadata), 10
axivity_set_session_id
 (axivity_get_session_id), 10
axivity_set_time (axivity_get_time), 11
axivity_unlock (axivity_is_locked), 11